Instructions for Use

For additional product information, visit www.vigilantbiosciences.com

Users must read this package insert in its entirety before using the product. Follow the instructions carefully when conducting the test. Failure to do so may cause inaccurate test results.

PRODUCT NAME AND INTENDED USE

The BeVigilant™ OraFusion System is to be used by qualified healthcare providers on patients 22 years of age or older to enhance the identification and screening of visible oral mucosal abnormalities. A qualitative screening result is generated by lateral flow immunoassays detecting two salivary biomarkers, p16 and EGFR, along with clinical risk factors. The BeVigilant™ OraFusion System is complimentary to and is intended to be used in combination with a traditional visual oral mucosal examination, as an aid in the assessment of oral cancer during routine examination. It is intended to assist in the decision regarding referral of the patient to a specialist. The BeVigilant™ OraFusion System is not automated.

PRINCIPLES OF OPERATION

The Ora-3D is a semi-quantitative immunoassay consisting of a single use sample The Ora-3D uses the platform of lateral flow immunoassay strips to detect p16 and EGFR in whole saliva, which are known to be associated with oral cancer. The strip includes a conjugate pad that is striped with antibodies, a test line that is striped with a capture antibody, and a control line that is striped with a 'catch-all' polyclonal antibody. The saliva specimen is added to the sample well of the cassette. When the saliva flows up the strip, the biomarkers bind with the antibody conjugates. This molecular complex binds to the antibody at the test line to give a positive signal, indicating presence of biomarkers in saliva. The control line develops a signal indicating that the strip performed properly.

The BeVigilant™ Reader, an electro-optical device, is used to read the Ora-3D test cassette and enter the clinical risk factors into the OraFusion Software. The Ora-3D test cassette is placed into the positioning tray on the BeVigilant™ Reader. The camera on the reader captures the image of the Ora-3D test cassette. The OraFusion Software processes the image and measures the critical line optical intensities. The absorbance units of the lines are used as an input to the algorithm, along with the clinical risk factors to generate a result.

Instructions for Use


MATERIALS PROVIDED IN THE BEVIGILANT™ READER TO ADMINISTER A TEST

Component Description (9000002)	Quantity
BeVigilant™ Reader	1
Power Adapter & Interchangeable Blade Kit	1
Instructions for Use (IFU)	1
Quick Reference Guide (QRG)	1

MATERIALS PROVIDED IN THE ORA-3D TO ADMINISTER A TEST

Component Description (9000006)	Quantity
Ora-3D Pouch	12
Oral Fluid Collection Device	12
Instructions for Use (IFU)	1
Quick Reference Guide (QRG)	1

Instructions for Use

ORA-3D & BEVIGILANT™ READER WARNINGS AND PRECAUTIONS

- The healthcare provider is required to complete training prior to the use of the BeVigilant™ OraFusion System. A Certification of Training will be provided to the healthcare provider. These required items can be found at www.vigilantbiosciences.com/training.
- Failure to follow the instructions provided may lead to inaccurate results.
- This test is not intended for diagnosis or as a standalone test.
- This test is intended for use by health-care providers in a health-care setting.
- The OraFusion System is intended for use on adults of age 22 years or older.
- This test has been developed for use with saliva only. The use of this test with any other specimen type may lead to inaccurate results.
- The use of a saliva specimen with blood, food, alcohol or any other substances present in the oral cavity within 1 hour may lead to false negative/positive results.
- Do not use this device if there is an infection or fungus in the mouth, blood present in the saliva sample, or if the patient has an active cancer diagnosis.
- Saliva sample should be used immediately after collection.
- Wear appropriate personal protective equipment and use standard office practices when handling and testing the patient specimen.
- Use the Ora-3D immediately after opening the pouch.
- Stop using immediately if an allergic reaction occurs.
- Insufficient saliva sample volume may lead to inaccurate results.
- This product has not been tested on pregnant women.
- This test is intended to be performed at room temperature (16°C to 25°C); do not use out of this range.
- The recommended storage temperature for the test kit is 4°C to 25°C.
- Use all test devices only once and dispose properly. Do not reuse any of the test devices.
- The BeVigilant™ Reader is not intended to be moved during a test.
- Do not use Ora-3D after the expiration date on the package. Expired devices should be disposed of according to local, state, and federal waste disposal requirements.
- Used Ora-3D are considered a potential biohazard and should be disposed of according to local, state, and federal waste disposal requirements.

Instructions for Use

- Do not position the equipment in a way that would make it difficult to reach the external power supply.
- Not properly plugging the charging cable into the BeVigilant™ Reader may cause the battery to lose charge and Reader to not power on.
- Ensure that the Reader is plugged into the appropriate power source.
- Do not expose ports to liquids; this can cause a short circuit and overheating.
- Maintain the environmental conditions for proper equipment operation.
- Modification of this device is not allowed.
- Ora-3D does not have backward compatibility.
- Use only with recommended consumables, accessories, or medical devices.
- Examine detachable power cord monthly for continued safe usage.
- The BeVigilant™ Reader may be reused until signs of material degradation occur. Do not use the BeVigilant™ Reader if the device shows signs of ageing, wear, fatigue, or any degradation as suggested by changes in its appearance that may affect performance.
- Install software updates immediately upon notification.
- Do not use Ora-3D if components are damaged or missing.
- Do not use the BeVigilant™ Reader if unresponsive or components are damaged or missing.
- Discontinue device use in the event of a malfunction.
- Label samples prior to testing to prevent mix-ups.
- Do not write over or damage the QR Code. If damage occurs, the cassette will return an invalid test result.
- Natural rubber latex was not used as a material in the manufacture of a medical product, its container and/or packaging.
- Avoid distractions while administering the test. The patient should remove any lipstick or lip treatment such as lip balm, ointment etc. prior to specimen collection for the administration of the test.

ORA-3D CONTRAINDICATIONS

There are no contraindications for the Ora-3D.

Instructions for Use

ORA-3D ENVIRONMENTAL CONDITIONS

The permissible environmental conditions of use for proper operation of the Ora-3D are:

Operational temperature:	16°C to 25°C
Operational Humidity:	10% to 85% relative humidity, non-condensing
Storage Temperature:	4°C to 25°C

BEVIGILANT™ READER PROCEDURE FOR INITIAL SET UP

1. Remove the BeVigilant™ Reader from the box.
2. Place the BeVigilant™ Reader on a flat open surface near a power outlet.
3. Remove the power adapter & Interchangeable Blade Kit from the box.
4. If needed, replace the interchangeable blade on the power adapter by pressing down on the spring-loaded tab, insert the blade adapter at 60° where the top is flat, and the bottom has a U-shape. Press the adapter down towards the power supply until it locks in place, a clicking sound will occur. If more information is needed, please visit: <http://www.globtek.com/pdf/Instructions-Interchangeable-Blades.pdf>
5. Plug into a power outlet.
6. Plug the power cord into the reader ensuring it is connected securely.
7. Press the power button, located on the back of the device. The power button will illuminate when device is on.
8. Ensure the BeVigilant™ Reader is powered on (if the screen does not light up after the power button is activated, check all connections and ensure the outlet being used is functioning properly).
9. Select the language to be used on the device.
10. Select the Time zone where the device is used.
11. Review the Privacy Policy and tap Accept. Tap CONTINUE to proceed.
12. Review the Terms of Use and tap Accept. Tap CONTINUE to proceed.
13. Connect the reader to a wireless internet network by selecting the network, entering the appropriate network password, and tap CONNECT.
14. Create an account for the practice by entering the Practice Name, Email Address, and Password (enter the password twice for confirmation) and tap CREATE ACCOUNT.
15. An automated message will be sent to the email address used to create the account. Check the email (if it is not in the inbox, check any spam filters) to get the Verification Code.

Instructions for Use

16. Enter the verification code in the reader. NOTE: Verification code is valid for a limited time.
17. If correct verification code is entered, the reader will prompt to create a 4-digit App PIN.
18. Create a 4-digit App PIN for the app to expedite future log in.
19. Re-enter the 4-Digit App PIN to confirm.
20. The reader is now ready to perform the Ora-3D. Set up is complete.
21. Please continue to the Vigilant website (<https://vigilantbiosciences.com/training/>) on a PC or mobile phone to complete training prior to use.

PLEASE NOTE: The USB is for Vigilant use in regard to internal testing and/or accessing data. Each Reader is password protected in order to prevent unauthorized access.

BEVIGILANT™ READER PROCEDURE FOR USE

Perform the following steps to initiate a new Ora-3D test:

1. Place the BeVigilant™ Reader on a flat open surface where the test will be performed. The reader should not be moved during testing. Movement may invalidate the test.
2. Tap the reader screen to ensure it is powered on. If the screen doesn't respond, press the power button on the back of the device for 1 second.
3. Enter the App PIN to log into the app.
4. Tap NEW TEST, then enter Year of Birth and Gender. Record the Test ID for traceability between the Ora-3D test results and the practice's Quality system for future reference.
5. Enter the patient Clinical Risk Factors by tapping the appropriate answer for each factor. (For additional information about Clinical Risk Factors, tap the ①).
6. Enter if the patient has fasted from food or fluid within 1 hour of the test. Enter if the patient has removed contaminants from lips. If the patient has not fasted the reader will show an error to instruct the patient to rinse with water and restart the test in one hour.
7. Prepare a saliva sample for testing by following the Specimen Collection Procedure for Use (next section). The BeVigilant™ OraFusion Software will act as a guide through the test steps.

Instructions for Use

SPECIMEN COLLECTION PROCEDURE FOR USE

Perform the following steps to prepare a patient sample for Ora-3D.

- Prior to collecting a sample, follow the BeVigilant™ Reader Procedure For Use.
- Gloves must be worn when performing the Ora-3D.

The OraFusion Software will guide the user through the following steps.

- Remove a Collection Device and a Test Cassette from the box. Inspect the packaging for damage and do not use items past the expiration date printed on the pouch. Only open the components as needed.


Important! Patients should not eat, drink, chew gum, smoke, gargle mouthwash, or put anything in or around the mouth for at least 1 hour before sample collection. Substances in the mouth may cause false test results.

IMPORTANT FOR SALIVA TRANSFER:


1. Place the sponge only on the top of the tongue. DO NOT place the sponge under the tongue.
2. Ensure to POOL saliva prior to placing the sponge in the mouth.
3. Maintain the tongue in one position while the sponge absorbs saliva.
4. DO NOT perform any sucking motion while the sponge is absorbing saliva.
5. DO NOT place your mouth over the plastic body of the collection device.

Instructions for Use

1. Instruct the patient to pool saliva in the mouth for a minimum of 30 seconds.
2. Remove the Push Cap and set aside.
3. Instruct the patient to place the entire Sponge in the mouth and saturate with saliva for a minimum of 90 seconds. Continue to pool and saturate the sponge in the mouth during 90 seconds.
4. After 90 seconds, check if the indicator has turned solid blue. If the indicator has not turned blue (indicating full saturation), repeat steps 1 through 3.
5. Retrieve the properly filled collection device from the patient. Place the Collection Device cap over the Sponge.
6. Keeping the device upright, push the cap down slowly and smoothly in one motion to filter the saliva into the sample dropper.
7. Saliva should be present in the bottom of the dropper tip. Ensure the minimum volume is up to the top of the purple dropper tip. If there is not enough, repeat steps 1 through 6.
8. Inspect the filtered saliva in the tube for blood or blood-like color. Select the option on the OraFusion Software to continue. If no blood/blood-like color is detected, then the sample is ready for testing. Proceed to the Ora-3D & BeVigilant™ Reader procedure for use.
9. Perform test within 1 hour of saliva collection.

Pool saliva in mouth for 30 seconds.


Remove Push Cap and set aside.


Put entire Sponge in mouth and saturate with saliva for 90 seconds. Do not chew or squeeze Sponge.

Indicator turns solid BLUE when saturated with saliva. Remove Sponge from mouth to check Indicator. If it is not solid Blue (fully activated), repeat steps 1 - 3.

Replace Push Cap over Sponge. Keep device upright.


Minimum volume of Saliva shall be to top of the dropper

Instructions for Use


ORA-3D & BEVIGILANT™ READER PROCEDURE FOR USE

1. Prior to administering a test, follow the BeVigilant™ Reader Procedure For Use and Specimen Collection Procedure For Use.
 2. Remove the Ora-3D Test cassette from the packaging and place on flat surface near the BeVigilant™ Reader and prepare to apply the filtered saliva to the cassette.
 3. Unscrew the dropper cap from the sample dropper nozzle.
- Caution:** Some saliva may flow without pressing the dropper tube. Ensure cleaning supplies are available if needed.
4. Add four drops of the sample into each well of the test cassette (eight drops in total) by squeezing the tube. The sample may pool in the well when first applied. Allow time for the saliva to saturate into the cassette before moving the cassette into the positioning tray.

Remove the dropper cap and set aside

Important! Do not exceed a 60 second timeframe between dropping the sample into the cassette and beginning the test.

Add a four (4) drops of the sample into each well of the test cassette by squeezing the tube.

5. Place the test cassette into the positioning tray on the right side of the BeVigilant™ Reader and slide the tray into the reader (do not lift Reader or positioning tray). The positioning tray and cassette have arrows indicating the direction of insertion.
6. Confirm all previous steps have been completed.
7. Tap START TEST on the BeVigilant™ Reader screen. If the elapsed time is greater than 1 hour since the lapse timer started, the sample has expired, and an error message will occur.
8. When the test is complete, the BeVigilant™ Reader will display a result. The result is based on the analysis of specific biomarkers and individual Clinical Risk Factors.

Instructions for Use

9. A result of "Monitor" suggests there is not an immediate need for specialist referral. A result of "Investigate Further" suggests further evaluation through a specialist referral.
10. Tap EMAIL RESULTS to email the test result. Tap BACK HOME to return to the home screen.
11. Remove the test cassette from the BeVigilant™ Reader and dispose of the cassette and Collection Device according to local, state, and federal waste disposal requirements.

If an error occurs, the screen will display an alert. A new test can be immediately administered using a new Ora-3D, a new collection device, and a fresh sample.

If you wish to view the test again, it can be displayed through the home screen by selecting "FIND TEST" and entering the Test ID previously recorded in the practice's Quality system.

INTERNAL QUALITY CONTROL

A control line is built into each Ora-3D, demonstrating the sample is properly flowing through the test. If the control line is not detected by the reader, the BeVigilant™ Reader will display an error message. A new test should be performed using a new collection device and fresh sample as explained above.

QUALITY CONTROL PROCEDURE

- Inspect the expiration date on the package for the Ora-3D to ensure the product is not expired.
- Examine detachable power cord monthly for continued safe usage.
- Ensure the most recent software is installed prior to use if the device has been unused for an extended period of time.

DISPOSAL

Used Ora-3D are considered a potential biohazard and should be disposed of according to local, state, and federal waste disposal requirements. Ensure that the BeVigilant™ Reader is handled in accordance with WEEE and Batteries Directive.

Instructions for Use

ANALYTICAL PERFORMANCE CHARACTERISTICS OF ORA-3D

The analytical performance results of the Ora-3D are shown in Table 1 below.

Table 1. Analytical Performance Characteristics

Analytical Performance Characteristic	p16 Result	EGFR Result
Analytical Sensitivity ¹	1.00	0.84
Trueness (Bias) ² PPA: Positive Percent Agreement NPA: Negative Percent Agreement	Accuracy. 92% PPA. 92% NPA. 91%	Accuracy. 88% PPA. 87% NPA. 91%
Linearity Range (Analytical Measuring Interval; AMI)	2.64 – 183 ng/mL	0.37 – 7.86 ng/mL
Measuring Range (Extended Measuring Interval; EMI)	N/A (limited to AMI)	10.82 – 150.82 ng/mL
Limit of Blank	0.18 ng/mL	0.15 ng/mL
Limit of Detection	2.64 ng/mL	0.37 ng/mL
Cut-Off ³	3.84 ng/mL	0.85 ng/mL
Precision (Repeatability) ANA: Average Negative Agreement APA: Average Positive Agreement	ANA. 97% APA. 98%	ANA. 94% APA. 97%
Analytical Performance Characteristic	OraFusion Exam	
Repeatability & Reproducibility ANA: Average Negative Agreement APA: Average Positive Agreement	ANA: 99% APA: 99%	
Analytical Specificity - Relevant Interferences	Refer to table 2.	
Limitations of Method	Saliva stability, sample delivery, chemical purity	
Available Reference Measurement Procedures	ELISA	

¹ The analytical sensitivity of a measuring system is the slope of the calibration curve in EP17-A2E. Calculations are based on standard curve calibrations using standard reference materials in whole human saliva. Calibration curves are established on a lot-to-lot basis.

² Based on measurement method comparison study vs. ELISA.

³ Ora-3D Assay cutoff values are based on clinical study results.

Instructions for Use

Table 2 Interference and Cross Reactivity Ora-3D

Compound	Test Dosage	Interference	Cross-Reactivity
<i>Pseudomonas aeruginosa</i>	1.00E+06 cfu/mL	Pass	Pass
<i>Streptococcus mutans</i>	1.00E+06 cfu/mL	Pass	Pass
<i>Eikenella corrodens</i>	1.00E+06 cfu/mL	Pass	Pass
<i>Candida albicans</i>	1.00E+06 cfu/mL	Pass	Pass
Herpes Simplex Virus	1.00E+05 pfu/mL	Pass	Pass
Epstein-Barr Virus	1.00E+05 pfu/mL	Pass	Pass
Cytomegalovirus	4.68E+03 pfu/mL	Pass	Pass
Human anti-mouse antibody	207 ng/mL	Fail	Pass
Topical Corticosteroid (Triamcinolone Acetonide Dental Paste)	5% v/v	Pass	Pass
Chlorhexidine Gluconate Oral Rinse (Peridex)	5% v/v	Pass	Pass
Human Serum Albumin	15 mg/mL	Pass	Fail
Immunoglobulin A	0.25 mg/mL	Pass	Pass
Lactoferrin	0.1 mg/mL	Pass	Pass
Coffee	5% v/v	Pass	Pass
Alcohol	5% v/v	Pass	Pass
Whole Human Blood	2.5% v/v	Pass	Pass
Toothpaste	0.5% w/v	Pass	Pass
Dental Adhesive / Mucoadhesive (Fixodent)	0.5% w/v	Pass	Pass
Mucin	2.5 mg/mL	Pass	Pass
Bilirubin	0.15 mg/mL	Pass	Pass
Alpha Amylase	250 U/mL	Pass	Pass
Alpha Amylase	250 U/mL	Pass	Pass
Anti-Inflammatory (Ibuprofen)	1060 U/mol/L	Pass	Pass
Topical Oral Pain Reliever (20% Benzocaine)	3 mg/mL / 0.3% w/v	Pass	Pass
Antiseptic Mouth Wash	5% v/v	Pass	Pass

CLINICAL PERFORMANCE CHARACTERISTICS OF ORA-3D+

The clinical performance characteristics of the Ora-3D++ are shown in Table 3 below.

Table 3. Clinical Performance Characteristics

Clinical Performance Characteristic	BeVigilant™ OraFusion System
Threshold Value	23.8%
Diagnostic Sensitivity	92.3%
Diagnostic Specificity	87.5%
Positive Predictive Value	64.9%
Negative Predictive Value	97.8%
Likelihood Ratio	[LR+] = 7.385 [LR-] = 0.088
Diagnostic Odds Ratio	84.0
Expected Values in Normal Populations	N/A
Expected Values in Affected Populations	N/A

+ These are ex-vivo, and known prevalence based whole human saliva samples.

++ These results are reflective of the BeVigilant™ OraFusion System.

Instructions for Use

GLOSSARY OF SYMBOLS

Standard/Source	Symbol	ISO/IEC symbol number	Title of Symbol	Description of Symbol Per Standard
ISO 15223-1:2021 5.1.1		ISO 7000-3082	Manufacturer	Indicates the medical device manufacturer
ISO 15223-1:2021 5.1.2		N/A	Authorized representative in the European Community/ European Union manufacturer	Indicates the authorized representative in the European Community/ European Union
ISO 15223-1:2021 5.1.3		ISO 7000-2497	Date of Manufacture	Indicates the date when the medical device was manufactured
ISO 15223-1:2021 5.1.4		ISO-7000-2607	Use-by-date	Indicates the date after which the medical device is not to be used.
ISO 15223-1:2021 5.1.5		ISO 7000-2492	Batch Code	Indicates the manufacturer's bath code so that the batch or lot can be identified.
ISO 15223-1:2021 5.1.6		ISO 7000-2493	Catalog number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
ISO 15223-1:2021 5.1.7		ISO 7000-2498	Serial Number	Indicates the manufacturer's serial number so that a specific medical device can be identified.
ISO 15223-1:2021 5.3.4		ISO 7000-0626	Keep Dry	Indicates a medical device that needs to be protected from moisture.

Instructions for Use

Standard/ Source	Symbol	ISO/IEC symbol number	Title of Symbol	Description of Symbol Per Standard
ISO 15223-1:2021 5.3.7		ISO 7000- 0632	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.
ISO 15223-1:2021 5.3.8		ISO 7000- 2620	Humidity limitation	Indicates the range of humidity to which the medical device can be safely exposed.
ISO 15223-1:2021 5.3.9		ISO 7000- 2621	Atmospheric pressure limitation	Indicates the range of atmospheric pressure to which the medical device can be safely exposed.
ISO 15223-1:2021 5.4.2		ISO 7000- 1051	Do not reuse	Indicates a medical device that is intended for one single use only.
ISO 15223-1:2021 5.4.3		ISO 7000- 1641	Consult instructions for use	Indicates that the user needs to consult the instructions for use.
ISO 15223-1:2021 5.5.1		N/A	In Vitro diagnostic medical device	Indicates a medical device that is intended to be used as an in vitro diagnostic medical device.
IVD Regulation 2017/746/EU		Annex I, Chapter III, 20.1h	Device for near-patient testing	Indication of near-patient testing
IVD Regulation 2017/746/EU		Annex I, Chapter III, 20.1h	Device not for self-testing	Explicit exclusion for assays not intended for self-testing or near-patient testing
IEC 60601- 1:2005 +AMD1:2012 +AMD2:2020		IEC 60417- 5032	Alternating current	Indicates that the equipment is suitable for alternating current only.

Instructions for Use

Standard/ Source	Symbol	ISO/IEC symbol number	Title of Symbol	Description of Symbol Per Standard
IEC 60601- 1:2005 +AMD1:2012 +AMD2:2020		IEC 60417- 5032	Alternating current	Indicates that the equipment is suitable for alternating current only.
IEC 60601- 1:2005 +AMD1:2012 +AMD2:2020		ISO 7010- W001	General warning sign	To signify a general warning
Waste from Electrical and Electronic Equipment (WEEE) Directive		IEC 60417- 6414	Waste Electrical and Electronic Equipment	Indicates that separate collection for waste electric and electronic equipment is required.
N/A		N/A	Illuminated Momentary Power Button	Indicates the power button on the back of the device. (Illuminated - Device is on / Light off - Device is off)
N/A		N/A	Saliva Drop Well indicator	Indicates where the user is to drop the saliva into the cassette.
N/A		N/A	CE Marking; European Conformity	Indicates that the manufacturer of the product complies with EU legislation and may be sold anywhere in the EEA (European Economic Area).
N/A		N/A	Prescription Use Only	Indicates prescription use only

Instructions for Use

Questions or Concerns?

If you have questions or concerns regarding these instructions for use, please contact your Sales Representative or Vigilant Biosciences at customerservice@vigilantbiosciences.com.

For additional product information, visit www.vigilantbiosciences.com

If any serious incident occurs in relation to the device, please report to Vigilant Biosciences, US Food & Drug Administration, and/or the competent authority of the Member State in which the user and/or the patient is established.

The Summary of Safety and Performance is available upon request.

REFERENCES

1. "The Science of Earlier: Improving early detection of oral and oropharyngeal cancer", White Paper, 2018, Vigilant Biosciences, Inc.
2. World Health Organization website. Accessed 2021.
3. Siegel RL, Miller KD, and Jemal A. Cancer statistics, 2019 CA Cancer J Clin 2019;69:7-34.

IVD For in *vitro* diagnostic use only.

PATIENT DATA PRIVACY

Vigilant Biosciences respects the privacy rights of individuals and is committed to handling and protecting personal information in compliance with the EU-U.S. and Swiss-U.S. Privacy Shield Frameworks. For our full Privacy Policy, visit our website at www.vigilantbiosciences.com. For complaints or concerns, contact our Privacy Administrator at privacy@vigilantbiosciences.com.

Instructions for Use

BeVigilant™ is a licensed trademark of Vigilant Biosciences, Inc.

©2024 by or under license to Vigilant Biosciences, Inc. All rights reserved. Subject to issued and pending international patent protection. All marks used herein are owned by Vigilant Biosciences, Inc. or their respective owners.

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Vigilant Biosciences, Inc.

1008 Ranch Road 620 South, Suite 204
Lakeway, TX 78734 USA
+1 800.979.4870
customerservice@VigilantBiosciences.com
www.vigilantbiosciences.com

Rx ONLY

Instructions for Use

Instructions for Use

Revision History

Document Revision	Description of Change	Effective Date
A	Initial Release	16-Apr-2021
B	Change in collection method	12-May-2021
C	Change in Authorized Representative	18-Apr-2022
D	Change in procedure to reflect product improvements	24-May-2022
E	Renaming of Device & Improvised Procedure Illustrations, Formatting correction, Warning addition	21-Mar-2023
F	Removal of duplicated warning, addition of Human serum albumin, coffee, mouthwash, toothpaste to a warning	28-Mar-2023
G	Update Age to Year of Birth, Addition of 'i' icon, Update to warning and precaution, Replace Wizard with OraFusion Software	26-APR-2023
H	Clarified the language in Intended Use, Principle of Operation, Addition of Warning related to training and patient age, removal of reader calibration check warning, Update to procedure to include training step, clarified patient fasting and contaminant step, clarified language in the check for blood step, update to results display, addition information about how to email results, Update to analytical and clinical performance characteristics. Updated Operational and Storage Temperature.	15-APR-2024
I	Clarified the language in Intended Use, Update to analytical and clinical performance characteristics.	22-NOV-2024

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BeVigilant™, Vigilant Biosciences, Inc.'in lisanslı bir ticari markasıdır.
Vigilant Biosciences GmbH